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

Dental System Control Panel

Home Page

- System Settings
 - System Settings
 - Design options
- DentalManager
 - General
 - Manufacturing Inbox
 - ERP and processing time
 - Order import
 - Working days
 - 3rd party applications
 - Milling machines
- Site Settings
 - This site
 - Dentists
 - Labs
 - Manufacturers
 - Operators
 - Countries
 - Manufacturing processes
- Tools
 - Subscription Management (Dongle)
 - Import/Export**
 - Download center
 - Training Center
 - 3Shape Communicate
- Basic elements
 - Materials
 - Colors
 - Margin line
 - Attachments
 - 2D Design overlays
 - CAD blocks
- Anatomy elements
 - ScanIt library
 - Anatomy and Pontic libraries
 - Smile libraries
 - Crowns
 - Crown Pontics
 - Inlays
 - Onlays/Veneers
 - Temporary Crowns
 - Temporary Pontics
 - Temporaries on prepared model
 - Tabletops
- Frame elements
 - Copings
 - Frame Pontics
 - Waxups
 - Primary telescopes
- Abutments
 - Top cap libraries
 - Implant systems
- Post and core
 - Post and core systems
 - Post and cores
- Bridge elements
 - Bars
 - Connectors
- Full dentures
 - Design settings
 - Gingivae
 - Smile libraries chart
- Removables
 - Wax profile strips
 - Retentions
 - Stippled waxes
- Digital model
 - Digital model design
 - Articulator interfaces
- Appliances
 - Orthodontics Control Panel

5. Seleccionar Import Materials.

Dental System Control Panel - Export/Import Materials

Home Back Save Help

Tools > Import/Export

Tasks

- Export materials
- Import materials

Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".

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
The screenshot shows the 'Dental System Control Panel - Export/Import Materials' interface. On the left, there are navigation buttons (Home, Back, Save, Help) and a 'Tools > Import/Export' section with 'Tasks' for 'Export materials' and 'Import materials'. The main area is titled 'Import materials' and contains instructions: 'To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".' Below the instructions is an 'Import materials' button. A 'History Of Import Materials' section is also visible.

Overlaid on this is a Windows Explorer window titled 'Abrir' (Open) showing the '3Shape' folder. The file list contains one entry:

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The 'Nombre:' field at the bottom of the Explorer window is empty, and the 'Abrir' (Open) button is highlighted.

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Dental System Control Panel - Implant system

Home Back Save Help Table mode Transform

Abutments > Implant systems

Implant system

- 3Shape Implant 4.3 mm
- B3i LDA (Certain/External)
- B3i LDA (Nobel Active)
- B3i LDA (Nobel Replace Select)
- B3i LDA (Straumann Bonelevel)
- DAS_3I_CTN_NP**
- DAS_E_C_1
- DAS_E_C_2

+ Add
+ Copy
X Delete
+ Move Up
+ Move Down
Recycle Bin
No items
+ Restore
Empty the Recycle Bin

Details
Total: 53 items

Edit Implant system

ID: 02049_DAS_3I_CTN_NP [Change ID](#)

Name: DAS_3I_CTN_NP

Version: V3

Abutment kits

- DAS_3I_CTN_NP_NR
- DAS_3I_CTN_NP_R

+ Add
+ Copy
X Delete
+ Edit
+ Preview

Implant models

- R3IINT34

+ Add
X Delete

Scan abutment models

- SBIND3NP

+ Add
X Delete

Screw models

- TPDH16L84

+ Add
X Delete

Abutment base models

- BASE_3I_CERTAIN_NP_NR
- BASE_3I_CERTAIN_NP_R

+ Add
X Delete

Abutment interface models

- IND3IHNP

+ Add
X Delete

Analog interface models