



# EXTRAORAL LIBRARY CODES

DYNAMIC ABUTMENT® SOLUTIONS LIBRARY



**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
BIOMET 3I	Certain	3,25/4	3,4	0,3	3,8/5/7/9	1,2	3,8
BIOMET 3I	Certain	4/5	4,1	0,3	4/5/7/9	1,2	4,0
BIOMET 3I	Certain	5	5	0,5	4,5		
BIOMET 3I	Osseotite	3,25	3,4	0,3	3,8/5/7/9	0,5	3,8
BIOMET 3I	Osseotite	4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
BIOMET 3I	Osseotite	5	5	0,5	4,5		
BIOMET 3I	Low profile	Universal	Universal	0,3	5,4/5/7/9		
AB	I2	3,5/3,75/4,2/4,5/5/6	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	I22	3,75/4,2	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	I5	3,5/3,75/4,2/4,5/5/6/7/8	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	I55	3,75/4,2/4,5/5/6/7/8	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	I10	4,2/5	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	I15	6/7/8	Standard	0,6	3,8/5/7/9	1,5	3,8
AB	Multi Unit D1P64	Universal	Universal	0,3	5,4/5/7/9		

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
2,0	3,8/5/7/9	3,0	3,8					DAS_E_0001
								DAS_E_0002
								DAS_E_0057
1,0	5/7/9							DAS_E_0003
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0058
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0025

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
ACE	External Hex	3,3	NP 3,5	0,3	3,8		
ACE	External Hex	3,75/4	RP 4,1	0,3	4/5/7/9	0,5	4/5/7/9
ACE	External Hex	4,75	WP 5	0,5	4,5		
ACE	Infinity Internal Hex	3,7/4,1	3,5	0,6	3,8/5/7/9	1,5	3,8
ACE	Infinity Internal Hex	4,7/5,1	4,5	0,4	4/5/7/9	1,0	4,0
ACE	Infinity TRICAM	3,5	3,5	0,5	3,8/5/7/9	1,2	3,8
ACE	Infinity TRICAM	4,3	4,3	0,3	4/5/7/9	1,2	4,0
ACE	Infinity TRICAM	5	5	0,3	4,0		
ACE	Infinity Octagon	3,3/4,1/4,8	RP 4,8	0,6	4,5		
ACE	Infinity Octagon	4,8	WP 6,5	0,6	4,5		
ACE	Multi Unit	Universal	Universal	0,3	5,4/5/7/9		
ADIN	Swell	3,3	3,45	0,6	3,8/5/7/9	1,5	3,8
ADIN	Swell	3,75/4,2	3,6	0,6	3,8/5/7/9	1,5	3,8
ADIN	Swell	5	4	0,6	3,8/5/7/9	1,5	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0058
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0026
								DAS_E_0027
								DAS_E_0028
								DAS_E_0037
								DAS_E_0096
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
ADIN	Swell	6	4,6	0,6	3,8/5/7/9	1,5	3,8
ADIN	TouregCloseFit	3	NP	1,5	4,5		
ADIN	TouregCloseFit	3,5	RP	1,5	3,8/5/7/9	2,0	3,8
ADIN	TouregCloseFit	4,3/5	WP	1,3	4/5/7/9	2,0	4,0
ADIN	TouregS/TouregOS	3,5	3,45	0,6	3,8/5/7/9	1,5	3,8
ADIN	TouregS/TouregOS	3,75/4,2	3,6	0,6	3,8/5/7/9	1,5	3,8
ADIN	TouregS/TouregOS	5	4	0,6	3,8/5/7/9	1,5	3,8
ADIN	TouregS/TouregOS	6	5	0,6	3,8/5/7/9	1,5	3,8
ADIN	Multi Unit TMA	Universal	Universal	0,3	5,4/5/7/9		
ALFA-GATE	Bioactive S-line/SLA	3,3/3,75/4,2/4,7/5,2/6	3,75	0,6	3,8/5/7/9	1,5	3,8
ALFA-GATE	Conical	3,5	NP	1,5	3,8/5/7/9	2,0	3,8
ALFA-GATE	Conical	4,3/5	RP	1,3	4/5/7/9	2,0	4,0
ALPHABIO	Conical Hex Connection (NICE)	3,2	Narrow	0,7	4,5	1,5	4,5
ALPHABIO	Conical Hex Connection (NEO)	3,2/3,5	Narrow	0,7	4,5	1,5	4,5

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0145
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022
		3,0	4,5					DAS_E_0136
		3,0	4,5					DAS_E_0136

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
ALPHABIO	Internal Hex Connection (IH) SPI	3,3/3,75/4,2/5/6	Universal	0,6	3,8/5/7/9	1,5	3,8
ALPHABIO	Internal Hex Connection (IH) ICE	3,7/3,75/4,2/4,65/5,3	Universal	0,6	3,8/5/7/9	1,5	3,8
ALPHABIO	Internal Hex Connection (IH) DFI	3,3/3,75/4,2/4,5	Universal	0,6	3,8/5/7/9	1,5	3,8
ALPHABIO	Internal Hex Connection (IH) ATID	3,3/3,75/4,2/5/6	Universal	0,6	3,8/5/7/9	1,5	3,8
ALPHABIO	Internal Hex Connection (IH) NEO	3,75/4,2/5	3,5	0,6	3,8/5/7/9	1,5	3,8
ALPHABIO	Conical Standard Connection (CS)	3,75/4,2/5	Standard	0,6	3,8	1,5	3,8
ANCLADEN	Anclalock	3,75/4,25/5	3,5	0,6	3,8/5/7/9	1,5	3,8
ANKYLOS	Ankylos	C/X	NP/RP/WP	1,0	3,8	1,0	3,8
ANTHOGYR	Axiom REG / PX	3,4	3,4 XNP	1,5	4,5		
ANTHOGYR	Axiom REG / PX	4	4 RP	1,5	4,0		
ANTHOGYR	Axiom REG / PX	4,6	4,6 WP	1,5	4,5		
ANTHOGYR	Anthofit HE	3,5/3,75/4	R (4,1)	0,3	4/5/7/9	0,5	4/5/7/9
ANTHOGYR	Ossfit	3,5/4,2	4,8	0,6	4,5		
ANTHOGYR	Ossfit	3,5/4,2	4,8	0,6	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
		3,0	3,8					DAS_E_0169
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0075
								DAS_E_0161
								DAS_E_0149
								DAS_E_0162
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0037
								DAS_E_0074

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
ANTHOGYR	Ossfit	5	6,5	0,6	4,5		
ANTHOGYR	Infinity Octagon	4,8	WP 6,5	0,6	4,5		
ANTHOGYR	Multi Unit	4,8	Universal	0,5	5,4		
ANTHOGYR	Axiom REG / PX	5,2	5,2	1,5	4,5		
ARDS	Classic	3,75/4,2/4,5	3,75	0,6	3,8/5/7/9	1,5	3,8
ARDS	Classic	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
ARDS	CIT	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
ARDS	Premium	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
ARDS	Smart	3,75/4,2/4,5	3,75	0,6	3,8/5/7/9	1,5	3,8
ASTRA	Yellow	3	Yellow	1,2	3,8		
ASTRA	Aqua	3,5/4	Aqua(Estrecha)	1,0	4/5/7/9	2,0	4,0
ASTRA	Lilac	4,5/5	Lilac (Ancha)	1,0	4/5/7/9	2,0	4,0
ASTRA	Cono 20°	Universal	Regular/Wide	0,6	3,8		
ASTRA	Evolution	3	3	1,0	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0096
								DAS_E_0096
								DAS_E_0163
								DAS_E_0162
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0109
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
								DAS_E_0066
								DAS_E_0090

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
ASTRA	Evolution	3,6	3,6	1,2	3,8	2,0	3,8
ASTRA	Evolution	4,2	4,2	1,5	4/5/7/9	2,0	4,0
ASTRA	Evolution	4,8	4,8	1,2	4,0	2,0	4,0
ASTRA	Evolution	5,4	5,4	1,0	4,5		
ASTRA	Evolution	Universal	Universal	0,5	3,8/5/7/9		
AVINENT	HE	3,3//3,5/4	3,5	0,3	3,8		
AVINENT	HE	3,3/3,8/4/4,2/4,8//4,5/5	4,1	0,3	4/5/7/9	0,5	4/5/7/9
AVINENT	HI	3,1//3,5/4	3,5	0,6	3,8/5/7/9	1,5	3,8
AVINENT	HI	3,3/3,8/4/4,2/4,8//4,5/5	4,1	0,6	3,8/5/7/9	1,5	3,8
AVINENT	Transepitelial	Universal	Regular Platform	0,3	5,4/5/7/9		
AVINENT	HE/EC	4,8	5,1	0,6	4,5		
BEGO	RS/RSX	3	3	0,6	3,8		
BEGO	S/RI/RS/RSX	3,25/3,75	3,25/3,75	0,7	3,8		
BEGO	S/RI/RS/RSX	4,1	4,1	0,5	4,0		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0	3,8	4,0	3,8					DAS_E_0006
								DAS_E_0007
								DAS_E_0091
								DAS_E_0092
								DAS_E_0008
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0025
								DAS_E_0061
								DAS_E_0049
								DAS_E_0050
								DAS_E_0051

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
BEGO	S/RI/RS/R SX	4,5	4,5	0,3	4,5		
BEGO	S/RI/RS/R SX	5,5	5,5	0,4	4,5		
BEGO	Hexágono Externo	2,7/2,9/3,1	Narrow Mini	0,3	3,8	0,5	3,8
BEGO	MultiPlus	Universal	Universal	0,5	5,4		
BIOCONCEPT	BC Tissue Level Standard	3,3/4,1/4,8	Regular	0,6	4,5		
BIOCONCEPT	BC Tissue Level Standard Plus	4,8	Regular	0,6	4,5		
BIOCONCEPT	BC Tissue Level Tapered Effect	4,8	Regular	0,6	4,5		
BIOCONCEPT	BC Bone Level	3,3	Narrow	1,3	3,8/5/7/9	2,0	3,8
BIOCONCEPT	BC Bone Level	4,1/4,8	Regular	1,1	4/5/7/9	2,0	4,0
BIOCONCEPT	BV Tapered Bone Level	3,5	Narrow	1,2	3,8/5/7/9	2,0	3,8/5/7/9
BIOCONCEPT	BV Tapered Bone Level	4/4,5/5	Regular	1,1	4/5/7/9	2,0	4,0
BIOGENESIS	3icon	3,3	Mini (Pink)	0,3	3,8		
BIOGENESIS	3icon	3,75/4/4,3/4,5	Regular (Blue)	0,3	4/5/7/9	0,5	4/5/7/9
BIOGENESIS	Aticon	3,5//4/4,5/5	Pink/Blue	1,0	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0052
								DAS_E_0081
1,0	3,8							DAS_E_0187
								DAS_E_0150
								DAS_E_0037
								DAS_E_0037
								DAS_E_0037
3,0		4,0		5,0				DAS_E_0033
3,0		4,0		5,0				DAS_E_0035
3,0	3,8	4,0	3,8					DAS_E_0029
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0	4,0	4,0	4,0					DAS_E_0005



**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
BIOGENESIS	Iticon	3,5/4,1/4,8	4,8	0,6	4,5		
BIOHORIZON	External	3,5	3,7 (Yellow)	0,3	3,8		
BIOHORIZON	Tapered Internal	3/3,4/3,8	3.0 (Grey)	1,8	3,8/5/7/9		
BIOHORIZON	Tapered Internal	3/3,8/4,6	3,5 (Yellow)	0,6	3,8/5/7/9	1,5	3,8
BIOHORIZON	Tapered Internal	4,6/5,8	4,5 (Green)	0,4	4/5/7/9	1,0	4,0
BIOHORIZON	Tapered Internal	5,8	5,7 (Blue)	0,4	4,5		
BIOHORIZON	Internal	3,5/4	3,5 (Yellow)	0,6	3,8/5/7/9	1,5	3,8
BIOHORIZON	Internal	4/5	4,5 (Green)	0,4	4/5/7/9	1,0	4,0
BIOHORIZON	Internal	5,8	5,7 (Blue)	0,4	4,5		
BIOHORIZON	Tapered Tissue Level	3/3,8	3,5 (Yellow)	0,6	3,8/5/7/9	1,5	3,8
BIOHORIZON	Tapered Tissue Level	4,6	4,5 (Green)	0,4	4/5/7/9	1,0	4,0
BIOHORIZON	Tapered Tissue Level	5,8	5,7 (Blue)	0,4	4,5		
BIOHORIZON	Multi Unit Angled	Universal	Universal	0,3	5,4/5/7/9		
BIOLOK	External Hex	3,45	3,45	0,3	3,8/5/7/9	0,5	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0037
								DAS_E_0023
								DAS_E_0102
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0080
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0080
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0080
								DAS_E_0025
1,0	5/7/9							DAS_E_0003

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
BIONER	Hikelt	3,8	3,95	0,6	3,8/5/7/9	1,5	3,8
BIONER	Hikelt	4,7	4,9	0,4	4/5/7/9	1,0	4,0
BIONER	Ikelt/Bikelt	3,3/3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
BIONER	Transepitelial A5M	Universal	Regular	0,3	5,4/5/7/9		
BIONER	TopDM	3,5	3,5	1,5	3,8/5/7/9	2,0	3,8
BIONER	TopDM	4	4	1,5	3,8/5/7/9	2,0	3,8
BIONER	TopDM	5	5	1,5	3,8/5/7/9	2,0	3,8
BIOTEC	SPR/CIM	3,3	3,3	0,6	3,8/5/7/9	1,5	3,8
BIOTEC	SPR/CIM	3,75	3,75	0,6	3,8/5/7/9	1,5	3,8
BIOTEC	SPR/SPTT/CIM	4,2	4,2	0,6	3,8/5/7/9	1,5	3,8
BIOTEC	SPR/SPTT/CIM	5	5	0,6	3,8/5/7/9	1,5	3,8
BIOTECH	Kontakt	3	Universal	1,5	4,5		
BIOTECH	Kontakt	3,6/4,2/4,8/5,4	Universal	1,0	4,0		
BREDDENT MEDICAL	Narrow Sky	3,5	NP 3,5	0,4	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0164
								DAS_E_0165
								DAS_E_0110

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
BREIDENT MEDICAL	Blue Sky	3,5/4/4,5/5,5	4	0,4	4,0		
BREIDENT MEDICAL	Blue Sky Classic	3,5/4/4,5	4	0,4	4,0		
BTI	External Connection	2,5/3/3,3/3,5/3,75	Tiny 3,5	0,3	3,8	0,5	3,8
BTI	External Connection	3,75/4/4,5/5	Universal 4,1	0,3	4/5/7/9	0,5	4/5/7/9
BTI	External Connection	4,5/5/5,5	Ancha 5,5	0,5	4,5		
BTI	Internal Connection	3,3/3,5/3,75/4/4,25/5/5,5	Universal 4,1	0,3	4,0		
BTI	Internal Connection	4,5/5/5,5/6/6,25	Ancha 5,5	0,7	4,5		
BTI	MultilM	Universal	Ancha 5,6	0,5	4,5	0,5	
BTK	Klassik / Konic EXT	3,25	3,4 EN	0,3	3,8/5/7/9	0,5	3,8
BTK	Klassik / Konic EXT	3,25PL/3,75/4	4,1 ER	0,3	4/5/7/9	0,5	4/5/7/9
BTK	Klassik / Konic INT	3,25/4	3,5 IR	0,6	3,8/5/7/9	1,5	3,8
B&W	Internal Connection	3,3/4	4	0,6	3,8		
B&W	External Connection	3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
CAMLOG	Conelog	3,8	3,8	1,0	4,0		

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0111
								DAS_E_0111
1,0								DAS_E_0009
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0060
								DAS_E_0010
								DAS_E_0059
								DAS_E_0177
1,0	5/7/9							DAS_E_0003
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0042
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0120

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
CAMLOG	Conelog	4,3	4,3	1,0	4,0		
CAMLOG	Screw Line 3,3	3,3	3,3	0,3	4,5		
CAMLOG	Screw Line 3,8	3,8	3,8	0,3	3,8		
CAMLOG	Screw Line 4,3	4,3	4,3	0,3	4,0		
CAMLOG	Screw Line 5,0	5	5	0,3	4,5		
CORTEX	Classix	3,3/3,8/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
CORTEX	Dynamix	3,3/3,8/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
CORTEX	Saturn	3,8/4,2	3,5	0,6	3,8/5/7/9	1,5	3,8
CORTEX	Conical Platform	3	NP	1,2	3,8		
CORTEX	Conical Platform	3,3/3,8/4,2	RP	1,0	4/5/7/9	2,0	4,0
CORTEX	Conical Platform	5/6	WP	1,0	4/5/7/9	2,0	4,0
CORTEX	Multi Unit	Universal	Universal	0,5	5,4/5/7/9		
CORTEX	Magix	3,3/3,8/4,2	RP	1,0	4/5/7/9	2,0	4,0
CONEXÃO SISTEMA DE PRÓTESE	Flash	3,5/4,3/5	Universal	1,5	3,8/5/7/9	2,0	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0121
								DAS_E_0087
								DAS_E_0011
								DAS_E_0012
								DAS_E_0088
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0109
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
								DAS_E_0025
3,0	4,0	4,0	4,0					DAS_E_0004
3,0		4,0		5,0				DAS_E_0021

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
CONEXÃO SISTEMA DE PRÓTESE	Torq	3,5/3,75/4	Universal	1,5	3,8/5/7/9	2,0	3,8
CONEXÃO SISTEMA DE PRÓTESE	Expand	3,75/4/5	Universal	1,5	3,8/5/7/9	2,0	3,8
DENTAURUM	Tiologic	3,3	Small	0,5	3,8		
DENTAURUM	Tiologic	3,7/4,2	Medium	0,5	4,0		
DENTAURUM	Tiologic	4,8/5,5	Large	0,5	4,5		
DENTAURUM	Tiiox	3,5	3,5	1,0	4,0		
DENTAURUM	Tiiox	4,5	4,5	1,0	4,0		
DENTEGRIS	SLS-Straight	4.5	4.5	0,4	4/5/7/9	1,0	4,0
DENTEGRIS	Sinus-Lift	4.5	4.5	0,4	4/5/7/9	1,0	4,0
DENTIUM	NR Line XNP	3,1	3,2	1,8	4,5		
DENTIUM	NR Line XNP	3,1	3,6	1,8	4,5		
DENTIUM	SuperLineImplantium	3,4	3,6	1,1	4/5/7/9	2,0	4,0
DENTIUM	SuperLineImplantium	3,8	4	1,1	4/5/7/9	2,0	4,0
DENTIUM	SuperLineImplantium	4,3	4,5	1,1	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0021
								DAS_E_0130
								DAS_E_0131
								DAS_E_0132
								DAS_E_0133
								DAS_E_0134
								DAS_E_0041
								DAS_E_0041
								DAS_E_0190
								DAS_E_0190
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
DENTIUM	SuperLineImplantium	4,8	5	1,1	4/5/7/9	2,0	4,0
DENTIUM	SuperLineImplantium	4,8	6	1,1	4/5/7/9	2,0	4,0
DENTIUM	SimpleLine II	3,8/4,3	4,8	0,6	4,5		
DENTIUM	SimpleLine II	3,8/4,3	4,8	0,6	4,5		
DENTIUM	SimpleLine II	4,3/4,8	6,5	0,6	4,5		
DENTIUM	NR LINE	3,6/4,3/5/6	3,6/4,3/5/6	0,8	3,8		
DENTIUM	NR Line Multi-Unit		5	0,3	5,4		
DENTIUM	Superline MU		4,5	0,3	5,1		
DENTIS	SClean	3,7	Mini	1,1	4/5/7/9	2,0	4,0
DENTIS	SClean	4,1/4,3	Regular	1,1	4/5/7/9	2,0	4,0
DENTIS	SClean	4,8	Wide	1,1	4/5/7/9	2,0	4,0
DIO IMPLANTS	External	3,3/3,8	Narrow 3,5	0,3	3,8		
DIO IMPLANTS	External	3,75/4/4,5	Regular 4,1	0,3	4/5/7/9	0,5	4/5/7/9
DIO IMPLANTS	External	5/5,3/5,5/6	Wide 5,1	0,6	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0037
								DAS_E_0074
								DAS_E_0096
								DAS_E_0191
								DAS_E_0192
								DAS_E_0193
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
DIO IMPLANTS	SM System	4,5/5/5,3	Regular/Wide	0,7	4,0		
DIO IMPLANTS	UF	3/3,3	Narrow	1,2	3,8	2,0	3,8
DIO IMPLANTS	UF	3,8/4/4,5/5/5,5	Regular	1,1	4/5/7/9	2,0	4,0
DIO IMPLANTS	Internal OCTA		4.8	0,6	4,5		
DSP BIOMEDICAL	External Hex	3,75/4/5//3,5/3,8/4,3	4,1	0,3	4/5/7/9	0,5	4/5/7/9
DENTAL TECH	Implogic	4,5	4,5 (Blue)	0,4	4/5/7/9	1,0	4,0
EASY IMPLANT	Multi Unit Conical Abutment	Universal	Universal	0,3	5,4/5/7/9		
EASY IMPLANT	Master C	3,5	3,5 (Ocean)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master C	4	4 (Ocean)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master C	4,5	4,5 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master C	5	5 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master S	3,3	3,3 (Ocean)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master S	3,75	3,75 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master S	4,25	4,25 (Lilas)	1,0	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0013
3,0		4,0		5,0				DAS_E_0014
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0074
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0041
								DAS_E_0025
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0	4,0	4,0	4,0					DAS_E_0004
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
EASY IMPLANT	Master S	4,75	4,75 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master L	3,3	3,3 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master L	3,75	3,75 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master L	4,25	4,25 (Lilas)	1,0	4/5/7/9	2,0	4,0
EASY IMPLANT	Master L	4,75	4,75 (Lilas)	1,0	4/5/7/9	2,0	4,0
ECKERMAN	Hexagon	3/3,5/4/4,5/5	4,1	0,3	4/5/7/9	0,5	4/5/7/9
ELITE MEDICAL	External Connection	3,75	Narrow Platform	0,3	3,8		
ELITE MEDICAL	External Connection	4	Regular Platform	0,3	4/5/7/9	0,5	4/5/7/9
ELITE MEDICAL	Conexión Externa	5	Wide	0,6	4,5		
EUROTEKNIKA	Aesthetica	4,1	4,8	0,6	4,5		
EUROTEKNIKA	Aesthetica	4,1	4,8	0,6	4,5		
EUROTEKNIKA	Aesthetica	4,8	6,5	0,6	4,5		
EUROTEKNIKA	Naturactis	3,5	3,4	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Naturactis	4	3,8	1,0	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0037
								DAS_E_0074
								DAS_E_0096
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004



**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
EUROTEKNIKA	Naturactis	4,5	4,3	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Naturactis	5	4,8	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Natea	3,6/4,1/4,8	Narrow	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Natea	3,6/4,1/4,8	Regular	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Natea	6	Wide	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Naturall	3,5	Narrow	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Naturall	4/4,5	Regular	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Naturall	5	Wide	1,0	4/5/7/9	2,0	4,0
EUROTEKNIKA	Uneva	3,6/4,1/4,8/6	4,1	0,3	4/5/7/9	0,5	4/5/7/9
EUROTEKNIKA	Multi Unit Tetra	Universal	Universal	0,3	5,4/5/7/9		
GALIMPLANT	Conexión interna	3,5	3,5	1,0	4/5/7/9	2,0	4,0
GALIMPLANT	Conexión interna	4	4	1,0	4/5/7/9	2,0	4,0
GALIMPLANT	Conexión interna	5	5	1,0	4/5/7/9	2,0	4,0
GALIMPLANT	Conexión externa	3,5/4	4	0,3	4/5/7/9	0,5	4/5/7/9

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0025
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
GALIMPLANT	Multi Unit Multiposicion recto	Universal	Universal	0,3	5,4/5/7/9		
GALIMPLANT	Multi Unit Multiposicion angulado	Universal	Universal	0,3	5,4/5/7/9		
GC TECH	AADVA Standard / Tapered Implants	3,3	NP	1,2	4,5	2,0	4,5
GC TECH	AADVA Standard / Tapered Implants	4	RP	1,2	3,8	2,0	3,8
GC TECH	AADVA Standard / Tapered Implants	5	WP	1,2	4,5		
GLOBAL D (TEKKA)	InKone Primo	3,5/4/4,5/5	5	1,0	4,5		
GLOBAL D (TEKKA)	InKone Universal	3,5/4/4,5/5	5	1,0	4,5		
GMI	Frontier	3,3/3,75/4,25	RP 3,3	0,6	3,8/5/7/9	1,5	3,8
GMI	Frontier	4,75/5,75	WP 4,3	0,4	4/5/7/9	1,0	4,0
GMI	Phoenix	3,3/3,75/4	Standard 4,1	0,3	4/5/7/9	0,5	4/5/7/9
GMI	Phoenix	5	Wide 5,1	0,6	4,5		
GMI	Pilar multiestético angulado	Universal	RP 4,8	0,3	5,4/5/7/9		
GT MEDICAL	Best Fit Internal Hex	3,7/4,1/4,3/4,8	Wide	1,0	4/5/7/9	2,0	4,0
GT MEDICAL	Best Fit External Hex	3,5	Narrow	0,3	3,8		

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0025
								DAS_E_0025
3,0	4,5							DAS_E_0196
3,0	3,8							DAS_E_0197
								DAS_E_0198
								DAS_E_0152
								DAS_E_0152
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0025
3,0	4,0	4,0	4,0					DAS_E_0005
								DAS_E_0023

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
GT MEDICAL	Best Fit External Hex	4,1	Regular	0,3	4/5/7/9	0,5	4/5/7/9
GT MEDICAL	Best Fit Hexágono Externo	5,1	Wide	0,6	4,5		
GT MEDICAL	Best Fit Internal Octagon	3,7/4,3/4,8	Regular	0,6	4,5		
GT MEDICAL	Best Fit Internal Octagon	3,7/4,3/4,8	Regular	0,6	4,5		
HANH IMPLANT	Tapered Implant	3,5/4,3	3,5 (Purple)	1,5	3,8/5/7/9	2,0	3,8
HANH IMPLANT	Tapered Implant	5	5 (Blue)	1,3	4/5/7/9	2,0	4,0
HANH IMPLANT	Tapered Implant	7	7 (Green)	1,4	4,5		
HANH IMPLANT	MultiUnit Abutment system	Universal	Universal	0,3	5,4/5/7/9		
HIOSSSEN	ETII SA / ETIII SA	3,5	Mini	1,2	3,8/5/7/9	2,0	3,8/5/7/9
HIOSSSEN	ETII SA / ETIII SA	4/4,5/5	Regular	1,1	4/5/7/9	2,0	4,0
HIOSSSEN	ETIII BA	3,5	Mini	1,2	3,8/5/7/9	2,0	3,8/5/7/9
HIOSSSEN	ETIII BA	4/4,5/5	Regular	1,1	4/5/7/9	2,0	4,0
HITEC	Tapered Self Thread	3,3/3,75	3,5	0,6	3,8/5/7/9	1,5	3,8
HITEC	Tapered Self Thread	4,2/5	4,5	0,4	4/5/7/9	1,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0037
								DAS_E_0074
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022
								DAS_E_0124
								DAS_E_0025
3,0	3,8	4,0	3,8					DAS_E_0029
3,0		4,0		5,0				DAS_E_0030
3,0	3,8	4,0	3,8					DAS_E_0029
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
HITEC	Logic Plus	3,5	3,7	0,6	3,8/5/7/9	1,5	3,8
HITEC	Logic Plus	4,3	3,9	0,6	3,8/5/7/9	1,5	3,8
IBS	Magic FC	4/4,5/5/5,5/6/6,5	3,8	1,1	4/5/7/9	2,0	4,0
IBS	N.R. Fix	3/3,5	3,8	1,1	4/5/7/9	2,0	4,0
IDO IMPLANTS	I Do Implants	3,8/4/4,5/5/5,5/6/7	Universal	1,1	4/5/7/9	2,0	4,0
IMPLANT DIRECT	Interactive	3,2/3,7	3	1,5	3,8/5/7/9	2,0	3,8
IMPLANT DIRECT	Interactive	4,3/5	3,4	1,3	4/5/7/9	2,0	4,0
IMPLANT DIRECT	Legacy	3,7	3,5	0,6	3,8/5/7/9	1,5	3,8
IMPLANT DIRECT	Legacy	4,7	4,5	0,4	4/5/7/9	1,0	4,0
IMPLANT DIRECT	ReplusReplantReactive	3,5/3,7/4,2	3,5	0,5	3,8/5/7/9	1,2	3,8
IMPLANT DIRECT	ReplusReplantReactive	4,3/4,7	4,3	0,3	4/5/7/9	1,2	4,0
IMPLANT DIRECT	ReplusReplantReactive	5/5,7	5	0,3	4,0		
IMPLANT DIRECT	SwishActive	3,3	3	1,5	3,8/5/7/9	2,0	3,8
IMPLANT DIRECT	SwishActive	4,1/4,8	3,4	1,3	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0026
								DAS_E_0027
								DAS_E_0028
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
IMPLANT DIRECT	SwishplantSwishplus	4,1/4,8	4,8	0,6	4,5		
IMPLANT DIRECT	SwishplantSwishplus	4,1/4,8	4,8	0,6	4,5		
IMPLANT DIRECT	Swishplant / Swishplus	4,8/5,7	6,5	0,6	4,5		
IMPLANT GENESIS	AKTIVE	3,5/3,75/4,2/5	Standard	0,6	3,8/5/7/9	1,5	3,8
IHDE DENTAL (IMBIODENT)	Bone Level Plus	3,3	3,3	1,3	3,8/5/7/9	2,0	3,8
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,1	4,1	1,1	4/5/7/9	2,0	4,0
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,8	4,8	1,1	4/5/7/9	2,0	4,0
INTRA LOCK	IntraHex	3,75/4	3,5	0,6	3,8/5/7/9	1,5	3,8
INTRA LOCK	IntraHex	4,75	4,5	0,4	4/5/7/9	1,0	4,0
INTRA LOCK	Unihex	4	Regular	0,3	4/5/7/9	0,5	4/5/7/9
INTRA LOCK	Unihex	4,75	Wide	0,3	4/5/7/9	0,5	4/5/7/9
JDENTALCARE	JDEvolution/JDEvolution Plus	3,7	3,7	0,6	3,8/5/7/9	1,5	3,8
JDENTALCARE	JDEvolution/JDEvolution Plus	4,3/5	4	0,6	3,8/5/7/9	1,5	3,8
JDENTALCARE	JDEvolution/JDEvolution Plus	6	5	0,6	3,8/5/7/9	1,5	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0037
								DAS_E_0074
								DAS_E_0096
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0033
3,0		4,0		5,0				DAS_E_0035
3,0		4,0		5,0				DAS_E_0035
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
JDENTALCARE	JD ICON	3,9	3,9	1,3	4/5/7/9	2,0	4,0
JDENTALCARE	JD ICON	4,3	4	1,3	4/5/7/9	2,0	4,0
JDENTALCARE	JD ICON	5	4,7	1,3	4/5/7/9	2,0	4,0
KEYSTONE	Internal TiLobe PrimaConnex	3,3/3,5	3,5	1,0	3,8/5/7/9		
KEYSTONE	Internal TiLobe PrimaConnex	4/4,1	4,1	1,0	4/5/7/9		
KEYSTONE	Internal TiLobe PrimaConnex	5	5	1,0	4/5/7/9		
KEYSTONE	Restore	3,75/4	RD 4,1	0,3	4/5/7/9	0,5	4/5/7/9
KLOCKNER	Essential Cone	3,5/4/4,5	4,5	0,6	4,0		
KLOCKNER	KL	3,5	Narrow	0,3	3,8		
KLOCKNER	KL	4,1	Regular	0,3	4/5/7/9	0,5	4/5/7/9
KLOCKNER	KL	5,1	Wide	0,6	4,5		
KLOCKNER	Vega	3,5	NV	1,2	3,8		
KLOCKNER	Vega	4/4,5	RV	1,2	4,0		
KLOCKNER	Multi Unit Permanent	4,2	Universal	0,5	5,4		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0022
3,0		4,0		5,0				DAS_E_0022
3,0		4,0		5,0				DAS_E_0022
								DAS_E_0044
								DAS_E_0045
								DAS_E_0046
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0054
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0082
								DAS_E_0083
								DAS_E_0173

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
LASAK	Bioniq	2,9	QN (Amarillo)	0,9	4,5		
LASAK	Bioniq	3,5/4/5	QR(Azul)	0,9	3,8		
LASAK	Multi Unit	Universal	Universal	0,3	5,4		
LEADER	Tixos External Hex	3,3/3,75	4,1	0,3	4/5/7/9	0,5	4/5/7/9
LEADER	Tixos Internal Hex	3,3	3,5	0,6	3,8/5/7/9	1,5	3,8
LEADER	Tixos Internal Hex	3,75	4	0,6	3,8/5/7/9	1,5	3,8
MEDENTIS (ICX)	Templant	3,75	3,75	1,1	4/5/7/9		
MEDENTIS (ICX)	Templant	4,1	4,1	1,1	4/5/7/9		
MEDENTIS (ICX)	Templant	4,8	4,8	1,1	4/5/7/9		
MEGAGEN	AnyOne External	3,5	Small 3,5	0,3	3,8		
MEGAGEN	AnyOne External	4	Regular 4,1	0,3	4/5/7/9	0,5	4/5/7/9
MEGAGEN	AnyOne External	4,5	Regular 4,5	0,3	4/5/7/9	0,5	4/5/7/9
MEGAGEN	AnyOne Internal	3,5/4/4,5/5/6/7	General	1,1	4/5/7/9	2,0	4,0
MEGAGEN	AnyRidge	3,5	Small	1,7	4/5/7/9	2,5	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0166
								DAS_E_0167
								DAS_E_0168
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0125
								DAS_E_0125
								DAS_E_0125
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0015

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
MEGAGEN	AnyRidge	4/4,5	Regular	1,7	4/5/7/9	2,5	4,0
MEGAGEN	AnyRidge	5/5,5	Wide	1,7	4/5/7/9	<b>2,5</b>	<b>4,0</b>
MEGAGEN	Mini	3/3,4	Mini	1,2	3,8	2,0	3,8
MEGAGEN	Cone Abutment	Universal	3,8	2,5	4,0		
MEGAGEN	Cone Abutment	Universal	4,8	0,6	4,5		
MEGAGEN	Multi Unit N Type	Universal	Universal	0,3	5,4/5/7/9		
MICRODENT	Ektos	3,7/4,2	3,5	0,6	3,8/5/7/9	1,5	3,8
MICRODENT	System	2,8/3,25	3,5	0,3	3,8/5/7/9	<b>0,5</b>	<b>3,8</b>
MICRODENT	Universal	2,8/3,25	3,5	0,3	3,8/5/7/9	0,5	3,8
MICRODENT	Universal	3,3/3,5/3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
MIS	C1	3,3	Narrow	0,8	3,8		
MIS	C1	3,75/4,2	Standard	0,7	4/5/7/9		
MIS	C1	5	Wide	1,2	4,0		
MIS	V3	3,9/4,3/5	Standard	0,7	4/5/7/9		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0015
								DAS_E_0015
3,0		4,0		5,0				DAS_E_0014
								DAS_E_0128
								DAS_E_0074
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
1,0	5/7/9							DAS_E_0003
1,0	5/7/9							DAS_E_0003
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0016
								DAS_E_0017
								DAS_E_0018
								DAS_E_0017



**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
MIS	Lance	3,75/4,2	Standard	0,3	4/5/7/9	0,5	4/5/7/9
MIS	M4	3,3	Narrow	0,8	3,8		
MIS	M4	3,75/4,2	Standard	0,6	3,8/5/7/9	1,5	3,8
MIS	M4	5/6	Wide	0,4	4/5/7/9	1,0	4,0
MIS	MultiUnit	Universal	General	0,6	5,1		
MIS	Seven	3,3	Narrow	0,8	3,8		
MIS	Seven	3,75/4,2	Standard	0,6	3,8/5/7/9	1,5	3,8
MIS	Seven	5/6	Wide	0,4	4/5/7/9	1,0	4,0
MOZOGRAU	MG Inhex	3,3	2,3 Mini	1,2	3,8		
MOZOGRAU	MG Inhex	3,75/4,25	2,8 Standard	1,0	4/5/7/9	2,0	4,0
MOZOGRAU	MG Inhex	5	3,8 Maxi	1,0	4/5/7/9	2,0	4,0
MOZOGRAU	MG Osseous	3,3	3,4 Mini	0,3	3,8/5/7/9	0,5	3,8
MOZOGRAU	MG Osseous	3,4/3,75/4,25	4,1 Standard	0,3	4/5/7/9	0,5	4/5/7/9
MOZOGRAU	MG Osseous	5	5 Maxi	0,6	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0019
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0020
								DAS_E_0019
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0109
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
1,0	5/7/9							DAS_E_0003
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
MPI	Privilege CM	3,5/4	Regular	1,0	4/5/7/9	2,0	4,0
MPI	Privilege CM	5	Wide	1,0	4/5/7/9	2,0	4,0
MPI	Excellence CM	3,5/4	Regular	1,0	4/5/7/9	2,0	4,0
MPI	Excellence CM	5	Wide	1,0	4/5/7/9	2,0	4,0
MPI	External Connection HE	3,3/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
NEOBIOTECH	EB External System	3,5	Narrow	0,3	3,8		
NEOBIOTECH	IS System	4	Regular 4	1,1	4/5/7/9	2,0	4,0
NEOBIOTECH	IS System	4,5	Regular 4,5	1,1	4/5/7/9	2,0	4,0
NEOBIOTECH	IS System	5	Wide 5	1,1	4/5/7/9	2,0	4,0
NEODENT	Cono Morse Mini Pilar	Universal	Regular	0,3	5,4/5/7/9		
NEODENT	Smart HE	3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
NEODENT	GM	3,5/3,75/4/4,3/5/6	Regular	1,2	4,0	2,5	4,0
NEOSS	ProActive	3,5 Green	ProActive	0,6	3,8/5/7/9		
NEOSS	ProActive	4 Yellow	ProActive	0,6	3,8/5/7/9		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0023
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0025
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,5	4							DAS_E_0186
								DAS_E_0047
								DAS_E_0047

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
NEOSS	ProActive	4,5 Blue	ProActive	0,6	4/5/7/9		
NEOSS	ProActive	5 Peach	ProActive	0,6	4/5/7/9		
NEOSS	ProActive	5,5 Lilac	ProActive	0,6	4/5/7/9		
NOBEL BIOCARE	Active	3	Mini 3.0	1,5	4,5		
NOBEL BIOCARE	Active	3,5	Narrow	1,5	3,8/5/7/9	2,0	3,8
NOBEL BIOCARE	Active	4,3/5	Regular	1,3	4/5/7/9	2,0	4,0
NOBEL BIOCARE	Active	5,5	Wide	1,4	4,5		
NOBEL BIOCARE	Branemark	3,3	Narrow	0,3	3,8		
NOBEL BIOCARE	Branemark	3,75/4	Regular	0,3	4/5/7/9	0,5	4/5/7/9
NOBEL BIOCARE	Branemark	5/6	Wide	0,6	4,5		
NOBEL BIOCARE	MultiUnit	4,8	Regular	0,3	5,4/5/7/9		
NOBEL BIOCARE	Replace	3,5	Narrow	0,5	3,8/5/7/9	1,2	3,8
NOBEL BIOCARE	Replace	4,3	Regular	0,3	4/5/7/9	1,2	4,0
NOBEL BIOCARE	Replace	5	Wide	0,3	4,0		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0048
								DAS_E_0048
								DAS_E_0048
								DAS_E_0159
3,0		4,0		5,0				DAS_E_0021
3,0		4,0		5,0				DAS_E_0022
								DAS_E_0124
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0025
								DAS_E_0026
								DAS_E_0027
								DAS_E_0028

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
NOBEL BIOCARE	Replace	6	Platform 6	0,3	4,5		
NOBEL BIOCARE	NobelSpeedy	3,3	Narrow	0,3	3,8		
NOBEL BIOCARE	NobelSpeedy	4/5	Regular	0,3	4/5/7/9	0,5	4/5/7/9
NOBEL BIOCARE	NobelSpeedy	5/6	Wide	0,6	4,5		
NORIS MEDICAL	TUFF	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
NORIS MEDICAL	TUFF TT	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
NORIS MEDICAL	ONIX	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
NORIS MEDICAL	CORTICAL	4/5/6.	3,75	0,6	3,8/5/7/9	1,5	3,8
NORIS MEDICAL	PTERYCORE	4,2	3,75	0,6	3,8/5/7/9	1,5	3,8
NORIS MEDICAL	PTERYFIT	4,2	3,75	0,6	3,8/5/7/9	1,5	3,8
NORMON	HI	3,75/4,25/4,75	3,5	0,6	3,8/5/7/9	1,5	3,8
NORMON	HE	3,25/3,75/4,25/4,75	4,1	0,3	4/5/7/9	0,5	4/5/7/9
NOVA IMPLANTS	PSI/PCI	3,3/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
OSSEOLIFE (ZIACOM)	OEX	3,75/4,25	RP 4,1	0,3	4/5/7/9	0,5	4/5/7/9

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0129
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
OSSTEM	TS	3,5	Mini	1,2	3,8/5/7/9	2,0	3,8/5/7/9
OSSTEM	TS	4/4,5/5/6/7	Regular	1,1	4/5/7/9	2,0	4,0
OSSTEM	US	3,3/3,5	Mini 3,5	0,3	3,8		
OSSTEM	US	3,75/4/4,5	Regular 4,1	0,3	4/5/7/9	0,5	4/5/7/9
OSSTEM	US	5/5,5	Wide 5,1	0,6	4,5		
OSTEOPLUS	SHE	3,45	3,45	0,3	3,8	0,5	3,8
OSTEOPLUS	SHE	3,75 / 4	4	0,3	4/5/7/9	0,5	4/5/7/9
OSTEOPLUS	SHI	3,3 / 3,75 / 4,2	3,5	0,6	3,8/5/7/9	1,5	3,8
OXY	K1 Line	3,5/4/4,5/5	Regular	1,7	4/5/7/9	2,5	4,0
OXY	PSK Line	3,5/4/4,5/5	Regular	1,7	4/5/7/9	2,5	4,0
OXY	MD Line KONE	4,75/4,25/5	Regular	1,7	4/5/7/9	2,5	4,0
OXY	MD Line Ext	3,75/4,25	Standard	0,3	4/5/7/9	0,5	4/5/7/9
PALTOP	ADVANCED CLASSIC	3,25	Narrow	0,5	4,5		
PALTOP	ADVANCED CLASSIC	3,75/4,2/5	Standard	0,6	3,8		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0	3,8	4,0	3,8					DAS_E_0029
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
1,0								DAS_E_0009
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0015
								DAS_E_0015
								DAS_E_0015
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0229
								DAS_E_0042

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
PALTOP	ADVANCED +	3,25	Narrow	0,5	4,5		
PALTOP	ADVANCED +	3,75/4,2/5	Standard	0,6	3,8		
PALTOP	DYNAMIC	3,25	Narrow	0,5	4,5		
PALTOP	DYNAMIC	3,75/4,2/5	Standard	0,6	3,8		
PALTOP	DIVA	3,75/4,2/5	Standard	0,6	3,8		
PALTOP	CONICAL ACTIVE	3,75/4,2/5	Standard	1,2	3,8/5/7/9	2,0	3,8/5/7/9
PALTOP	Universal Multi-Unit		Universal	0,4	5,4/5/7/9		
PHIBO	BTN Serie 3	3,6	4	0,3	4/5/7/9	0,5	4/5/7/9
PHIBO	BTN Serie 4	4,2	4	0,3	4/5/7/9	0,5	4/5/7/9
PHIBO	TSH Serie 3	3,6	4	0,3	4/5/7/9	0,5	4/5/7/9
PHIBO	TSH Serie 4	4,2	4	0,3	4/5/7/9	0,5	4/5/7/9
PROCLINIC	Cilindrico Externo/Cónico Externo	3,75/4,25//3,5/4	4,1 Estandar	0,3	4/5/7/9	0,5	4/5/7/9
PROCLINIC	Cilindrico Interno/Cónico Interno	3,3/3,75/4,25/5//3,5/4/5	3,5	0,6	3,8/5/7/9	1,5	3,8
PROCLINIC	SP Octa	3,3/4,1/4,8	4,8	0,6	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0229
								DAS_E_0042
								DAS_E_0229
								DAS_E_0042
								DAS_E_0042
3,0	3,8	4,0	3,8					DAS_E_0029
								DAS_E_0181
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0037

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
PROCLINIC	SP Octa	3,3/4,1/4,8	4,8	0,6	4,5		
PROCLINIC	SP Octa	4,8	6,5	0,6	4,5		
PROCLINIC	Aqua CM	3,5/4/5	2,82	1,0	4/5/7/9	2,0	4,0
RADHEX	PHE	3,5	3,5	0,3	3,8		
RADHEX	PHE	4/4,5/5	4,1	0,3	4/5/7/9	0,5	4/5/7/9
RADHEX	PHI	3,75	3,5	0,6	3,8/5/7/9	1,5	3,8
RADHEX	PHI	4,5/5	4,5	0,4	4/5/7/9	1,0	4,0
SEWON MEDIX	IH2 SLA SYSTEM	3,5	Mini	1,2	3,8/5/7/9	2,0	3,8/5/7/9
SEWON MEDIX	IH2 RBM SYSTEM	3,5	Mini	1,2	3,8/5/7/9	2,0	3,8/5/7/9
SEWON MEDIX	IH2 SLA SYSTEM	3,5/4/4,5/5	Regular	1,1	4/5/7/9	2,0	4,0
SEWON MEDIX	IH2 RBM SYSTEM	3,5/4/4,5/5	Regular	1,1	4/5/7/9	2,0	4,0
SEWON MEDIX	IH SYSTEM	3,5/4/4,5/5	Universal	0,3	5,4/5/7/9		
SGS	P1	3,2/3,75/4,2/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8
SGS	P7	3,2/3,75/4,2/4,5/5/6	3,75	0,6	3,8/5/7/9	1,5	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0074
								DAS_E_0096
3,0	4,0	4,0	4,0					DAS_E_0004
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
3,0	3,8	4,0	3,8					DAS_E_0029
3,0	3,8	4,0	3,8					DAS_E_0029
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
3,0		4,0		5,0				DAS_E_0040

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
SIC INVENT	HEXAGONAL SYSTEM SICace	3,4/4	3,3	0,3	4,5/5/7/9		
SIC INVENT	HEXAGONAL SYSTEM SICace	4,5/5	4,2	0,3	3,8		
SIGNO VINCES	Duo	4,6	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SIGNO VINCES	Inttegra	3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SIGNO VINCES	Compact	4,5	CM3,8	1,0	4/5/7/9	2,0	4,0
SIGNO VINCES	Duocon	3,8	CM3,8	1,0	4/5/7/9	2,0	4,0
SIGNO VINCES	Duocon	4,6/5,5	CM4,6	1,0	4/5/7/9	2,0	4,0
SIGNO VINCES	Infra	3,3/3,8/4,6	CM	1,0	4/5/7/9	2,0	4,0
SIN IMPLANTS	HI SW	3,8	3,8	0,7	4/5/7/9	0,7	4,0
SIN IMPLANTS	Tryon	3,25/3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SIN IMPLANTS	Tryon CO	4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SIN IMPLANTS	Revolution	3,25/3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SIN IMPLANTS	Stylus	4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SOUTHERN IMPLANTS	TriNex	3,5	3,5	0,5	3,8/5/7/9	1,2	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0170
								DAS_E_0171
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0004
								DAS_E_0039
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0026



**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
SOUTHERN IMPLANTS	TriNex	4,3	4,3	0,3	4/5/7/9	1,2	4,0
SOUTHERN IMPLANTS	TriNex	5	5	0,3	4,0		
SOUTHERN IMPLANTS	TriNex	6	6	0,3	4,5		
SOUTHERN IMPLANTS	IT Connection	3,3/4/4,1/4,9/5	4,8	0,6	4,5		
SOUTHERN IMPLANTS	IT Connection	4,9/5/6	6,5	0,6	4,5		
SOUTHERN IMPLANTS	External Hex	3,25	3,4	0,3	3,8/5/7/9	0,5	3,8
SOUTHERN IMPLANTS	External Hex	3,75/4	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SOUTHERN IMPLANTS	Internal Hex	3,75/4,2/5	Universal	0,6	3,8/5/7/9	1,5	3,8
SOUTHERN IMPLANTS	Compact Conical MultiUnit	4,8	4,8	0,3	5,4/5/7/9		
SOUTHERN IMPLANTS	Deep Conical	3	2,45	1,2	3,8		
SOUTHERN IMPLANTS	Deep Conical	3,5/4	2,95/3,1	1,0	4/5/7/9	2,0	4,0
SOUTHERN IMPLANTS	Deep Conical	5	4,1	1,0	4/5/7/9	2,0	4,0
STERIOSS	HexLoc	3,25	3,3	0,3	3,8		
STERNGOLD	Stern Ex	3,75/4/5	4,1	0,3	4/5/7/9	0,5	4/5/7/9

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0027
								DAS_E_0028
								DAS_E_0129
								DAS_E_0037
								DAS_E_0096
1,0	5/7/9							DAS_E_0003
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0025
								DAS_E_0109
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
SYBRON	Endopore	4,1	4,1	0,3	4/5/7/9	0,5	4/5/7/9
SWEDEN MARTINA	Premium khono	3,3	3,3	0,5	3,8		
SWEDEN MARTINA	Premium khono	3,8	3,8	0,7	4,0		
SWEDEN MARTINA	Premium khono	4,25	4,25	0,9	4,0		
SWEDEN MARTINA	Outlink 2	3,3	3,3	0,6	3,8		
SWEDEN MARTINA	Outlink 2	3,75/4,1	4,1	0,6	4,0		
SWEDEN MARTINA	Prama	3,8/4,25/5	3,4	0,5	3,8		
SWEDEN MARTINA	Shelta	3,8	3,8	0,7	4,0		
SWEDEN MARTINA	Shelta	4,25	4,25	0,7	4,0		
SWEDEN MARTINA	Shelta	5	5	0,7	4,0		
STRAUMANN	Standard Plus	3,3	NNC 3,5	0,5	4,5		
STRAUMANN	Bone Level	3,3	NC 3,3	1,3	3,8/5/7/9	2,0	3,8
STRAUMANN	Bone Level	4,1/4,8	RC4,1/4,8	1,1	4/5/7/9	2,0	4,0
STRAUMANN	BLX	3,5/3,75/4/4,5	RB (Regular Base)	1,5	4,5/5/7/9		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0031
								DAS_E_0032
								DAS_E_0065
								DAS_E_0063
								DAS_E_0064
								DAS_E_0031
								DAS_E_0032
								DAS_E_0032
								DAS_E_0032
								DAS_E_0160
3,0		4,0		5,0				DAS_E_0033
3,0		4,0		5,0				DAS_E_0035
3,0	5/7/9							DAS_E_0207

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
STRAUMANN	BLX	5/5,5/6,5	WB (Wide Base)	1,5	4,5		
STRAUMANN	Internal Octagon	3,3/4,1/4,8	Regular 4,8	0,6	4,5		
STRAUMANN	Tissue Level	4,8	Wide 6,5	0,6	4,5		
STRAUMANN	Synocta	Universal	Regular 4,8	0,6	4,5		
STRAUMANN	Bone Level Tapered SC	2,9	SC 2,9	1,0	4,5		
STRAUMANN	Synocta	Universal	Wide 6,5	0,6	4,5		
STRAUMANN	Screw Retained	Universal	NC/RC	0,5	5,2		
STRAUMANN	Tissue Level	3,3	3,5	0,5	4,5		
TBR	Hexágono Externo	3,5	Narrow	0,3	3,8		
TBR	Hexágono Externo	5	Wide	0,5	4,5		
	Octógono Interno	4	4	0,5	3,8		
TITANIUMFIX	Bfix	3,5/4	Regular	1,0	4/5/7/9	2,0	4,0
TITANIUMFIX	Bfix	4,5/5	Larga	1,0	4/5/7/9	2,0	4,0
TREE-OSS	Simple	3,3/3,75/5	3,75	0,6	3,8/5/7/9	1,5	3,8

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0208
								DAS_E_0037
								DAS_E_0096
								DAS_E_0074
								DAS_E_0135
								DAS_E_0137
								DAS_E_0101
								DAS_E_0160
								DAS_E_0023
								DAS_E_0058
								DAS_E_0267
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0		4,0		5,0				DAS_E_0040

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
TREE-OSS	Hexagono Externo	3,3	3,5 Rosa	0,3	3,8		
TREE-OSS	Hexagono Externo	3,75/4	4,1 Amarillo	0,3	4/5/7/9	0,5	4/5/7/9
TREE-OSS	Hexagono Externo	5	5,1 Azul	0,6	4,5		
TREE-OSS	Tri-Channel Interno	3,5	3,5 Rosa	0,5	3,8/5/7/9	1,2	3,8
TREE-OSS	Tri-Channel Interno	4,3	4,3 Amarillo	0,3	4/5/7/9	1,2	4,0
TREE-OSS	Tri-Channel Interno	5	5,0 Azul	0,3	4,0		
TREE-OSS	MULTI UNIT	Universal	Universal	0,3	5,4/5/7/9		
TRI DENTAL IMPLANTS	TriVent	3,75/4,1/4,7	3,5	0,6	3,8/5/7/9	1,5	3,8
TRINON	Q2	3,5/3,75/4,5	4	0,3	4/5/7/9	0,5	4/5/7/9
TRINON	QK	4	4,8	0,6	4,5		
TRINON	QK	4	4,8	0,6	4,5		
UFIT	Gt2	3,5	Mini	1,0	4/5/7/9	2,0	4,0
UFIT	Gt2	4/4,5	Regular	1,0	4/5/7/9	2,0	4,0
UFIT	Gt2	5	Wide	1,0	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
								DAS_E_0023
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0061
								DAS_E_0026
								DAS_E_0027
								DAS_E_0028
								DAS_E_0025
3,0		4,0		5,0				DAS_E_0040
1,0	5/7/9	2,0	5/7/9	3,0	5/7/9	4,0	5/7/9	DAS_E_0024
								DAS_E_0037
								DAS_E_0074
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0005

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
UFIT	Gt2	5,5/6/6,5/7	Ultrawide	1,0	4/5/7/9	2,0	4,0
UFIT	Gt2	3,5	Mini	1,0	4/5/7/9	2,0	4,0
UFIT	Nt2	3,5	Mini	1,0	4/5/7/9	2,0	4,0
UFIT	Nt2	4/4,5	Regular	1,0	4/5/7/9	2,0	4,0
UFIT	Nt2	5	Wide	1,0	4/5/7/9	2,0	4,0
UFIT	Nt2	5,5/6/6,5/7	Ultrawide	1,0	4/5/7/9	2,0	4,0
VULKAN IMPLANTS	INHex	3,3/3,75/4,2/5	3,75	0,6	3,8/5/7/9	1,5	3,8
XIVE	S	3,4	3,4	0,7	3,8/5/7/9	0,7	3,8
XIVE	S	3,8	3,8	0,7	4/5/7/9	0,7	4,0
XIVE	S	4,5	4,5	0,3	4,0		
XIVE	XIVE	3	3	1,0	4,5		
XIVE	XIVE	5,5	5,5	1,0	4,5		
YES IMPLANT	SSYSTEM	3,3/3,5	Narrow	1,1	4/5/7/9	2,0	4,0
YES IMPLANT	SSYSTEM	4/4,5	Regular	1,1	4/5/7/9	2,0	4,0

**NOTE:**

DAS = Dynamic Abutment Solutions  
E = Extraoral  
NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
GH = Gingiva Height  
CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0004
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0005
3,0	4,0	4,0	4,0					DAS_E_0005
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0038
								DAS_E_0039
								DAS_E_0085
								DAS_E_0084
								DAS_E_0086
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



BRAND	MODEL	Ø IMPLANT	PLATFORM	GH A	CH A	GH B	CH B
YES IMPLANT	SSYSTEM	5/5,5	Wide	1,1	4/5/7/9	2,0	4,0
YES IMPLANT	SSYSTEM	6/6,5/7/8	BWide	1,1	4/5/7/9	2,0	4,0
ZIMMER	SwissPlus	3,7/4,1/4,8	4,8	0,6	4,5		
ZIMMER	SwissPlus	3,7/4,1/4,8	4,8	0,6	4,5		
ZIMMER	Zimmer Type (Internal Hexagon)	3,7/4,1	3,5	0,6	3,8/5/7/9	1,5	3,8
ZIMMER	Zimmer Type (Internal Hexagon)	4,7	4,5	0,4	4/5/7/9	1,5	4,0
ZIMMER	ScrewVent	3,7/4,1	3,5	0,6	3,8/5/7/9	1,5	3,8
ZIMMER	ScrewVent	4,7	4,5	0,4	4/5/7/9	1,5	4,0
ZIMMER	ScrewVent	6	5,7	0,4	5,7		
ZIMMER	ScrewVent	6	5,7	0,4	5,7		
ZIMMER	Tapered	Universal	Universal	0,3	3,8		
ZIMMER	Eztetic	3,1	2,9	1,5	4,5		

**NOTE:**

DAS = Dynamic Abutment Solutions  
 E = Extraoral  
 NR = Non-rotator (Engaging)

R = Rotator (Non-Engaging)  
 GH = Gingiva Height  
 CH = Cement Height



GH C	CH C	GH D	CH D	GH E	CH E	GH F	CH F	LIBRARY CODE
3,0		4,0		5,0				DAS_E_0030
3,0		4,0		5,0				DAS_E_0030
								DAS_E_0037
								DAS_E_0074
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
3,0		4,0		5,0				DAS_E_0040
								DAS_E_0041
								DAS_E_0056
								DAS_E_0080
								DAS_E_0205
								DAS_E_0178